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Method 1

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), Nonselective Sheep Blood Agar, Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 2

Media	Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 72 hours

Method 3

Media	Chocolate Agar
Temperature	35°C
Atmosphere	5 to 7% Carbon Dioxide
Growth Time	24 to 48 hours

Method 4

Media	Anaerobic Blood Agar
Temperature	35°C
Atmosphere	Anaerobic
Growth Time	48 to 72 hours

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Method 6

Media	Chocolate Agar
Temperature	35°C
Atmosphere	Microaerophilic
Growth Time	48 to 72 hours

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Method 7

Media	Lowenstein Jensen Agar or Middlebrook Agar
Temperature	35°C
Atmosphere	5 to 7% CO ₂ or aerobic atmosphere.
Growth Time	2 to 30 days

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Method 8

Media	Buffered Charcoal Yeast Extract Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	3 to 5 days

Method 9

Media	V Agar or Chocolate Agar
Temperature	35°C
Atmosphere	5 to 7% Carbon Dioxide
Growth Time	48 hours

Method 10

Media	Tryptic Soy Agar (Soybean Casein Digest Agar) or Marine Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

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Method 11

Media	Phase 1: MRS (Man, Rogosa, Sharpe) Broth Phase 2: Columbia CNA with Sheep Blood or Tryptic Soy Agar with Sheep Blood
Temperature	Phase 1: 35°C Phase 2: 35°C
Atmosphere	Phase 1: Aerobic Phase 2: 5 to 7% Carbon Dioxide
Growth Time	Phase 1: 48 hours Phase 2: 48 hours

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Method 12

Media	Potato Dextrose Agar
Temperature	55°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 13

Media	Lowenstein Jensen Agar or Middlebrook Agar
Temperature	35°C
Atmosphere	5 to 7% CO ₂ or Aerobic
Growth Time	May require up to one month of incubation.

Method 14

Media	Leeming Notman Agar
Temperature	30°C
Atmosphere	Aerobic
Growth Time	72 hours

Method 15

Media	Chocolate agar, Sheep Blood Agar, Tryptic Soy Agar, and Bordet Gengou Agar with 15% Defibrinated Sheep Blood
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

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Method 16

Media	Chocolate agar or Bordet Gengou Agar with 15% Defibrinated Sheep Blood
Temperature	35°C
Atmosphere	Aerobic
Growth Time	2 days to 1 week

Method 17

Media	Phase 1: ISF (modified Infant Soy Formula) Broth Phase 2: Sulfite Agar
Temperature	Phase 1: 55°C Phase 2: 55°C
Atmosphere	Phase 1: Anaerobic Phase 2: Anaerobic
Growth Time	Phase 1: 48 hours Phase 2: 48 hours to 7 days

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Method 18

Media	Middlebrook 7H11 Agar
Temperature	30°C
Atmosphere	5 to 7% CO ₂
Growth Time	3 to 4 weeks

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Method 19

Media	Sheep Blood Agar supplemented with Pyridoxal
Temperature	35°C
Atmosphere	5 to 7% CO ₂
Growth Time	24 to 48 hours

Method 20

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), nonselective Sheep Blood Agar, Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	6 days

Method 21

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), nonselective Sheep Blood Agar, Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 22

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), nonselective Sheep Blood Agar, Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	30°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 23

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), Nonselective Sheep Blood Agar or Nutrient Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 24

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), Nonselective Sheep Blood Agar, Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	55°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 25

Media	Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	5 days

Method 26

Media	Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	5 days

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Method 27

Media	Tryptic Soy Agar (Soybean Casein Digest Agar) or Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 28

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), Nonselective Sheep Blood Agar, or Standard Methods Agar (Plate Count Agar)
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 29

Media	Nonselective Sheep Blood Agar, Standard Methods Agar (Plate Count Agar), or Nutrient Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 30

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), Nonselective Sheep Blood Agar, Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

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Method 31

Media	Nonselective Sheep Blood Agar
Temperature	30°C
Atmosphere	Aerobic
Growth Time	24 to 72 hours

Method 32

Media	Nonselective Sheep Blood Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	24 to 72 hours

Method 33

Media	Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	5 to 10% CO ₂
Growth Time	24 to 72 hours

Method 34

Media	Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 72 hours

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Method 35

Media	Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	Microaerophilic
Growth Time	24 to 72 hours

Method 36

Media	Chocolate Agar
Temperature	25°C
Atmosphere	5 to 7% Carbon Dioxide
Growth Time	3 to 4 days

Method 37

Media	Chocolate Agar
Temperature	35°C
Atmosphere	5 to 7% Carbon Dioxide
Growth Time	24 to 48 hours

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Method 38

Media	Anaerobic Blood Agar
Temperature	35°C
Atmosphere	Anaerobic
Growth Time	2 to 5 days

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Method 42

Media	Anaerobic Blood Agar
Temperature	35°C
Atmosphere	Anaerobic
Growth Time	3 to 4 days

Method 43

Media	Anaerobic Blood Agar
Temperature	35°C
Atmosphere	Anaerobic
Growth Time	5 to 7 days

Method 44

Media	Anaerobic Blood Agar
Temperature	35°C
Atmosphere	Anaerobic
Growth Time	3 to 5 days

Method 45

Media	Potato Dextrose Agar
Temperature	45°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 46

Media	Sabouraud Dextrose Emmons Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	2 to 7 days

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Method 47

Media	Sabouraud Dextrose Emmons Agar or Malt Extract Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	2 to 7 days

Method 48

Media	Potato Dextrose Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	7 days

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Method 50

Media	Sabouraud Dextrose Emmons Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	2 to 7 days

Method 51

Media	Sabouraud Dextrose Emmons Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	7 to 14 days

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Method 54

Media	See notes below for rehydration instructions; Marine Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

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Method 55

Media	See note below for important directions. Phase 1: MRS (Man, Rogosa, Sharpe) Broth Phase 2: MRS Agar
Temperature	Phase 1: 25°C Phase 2: 25°C
Atmosphere	Phase 1: Aerobic Phase 2: 5 to 7% Carbon Dioxide
Growth Time	Phase 1: 48 to 72 hours Phase 2: 72 to 96 hours.

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Method 56

Media	See note below for important directions. Phase 1: Cooked Meat Medium Phase 2: Anaerobic Blood Agar
Temperature	Phase 1: 45°C Phase 2: 45°C
Atmosphere	Phase 1: Aerobic Phase 2: Anaerobic
Growth Time	Phase 1: 72 hours Phase 2: 3 to 5 days

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Method 57

Media	Bordet Gengou Agar with 15% Defibrinated Sheep Blood
Temperature	35°C
Atmosphere	Aerobic
Growth Time	2 days to 1 week

Method 58

Media	See note below for important directions. Phase 1: Mycoplasma Broth Phase 2: Mycoplasma Agar
Temperature	Phase 1: 35°C Phase 2: 35°C
Atmosphere	Phase 1: Aerobic Phase 2: 5 to 7% CO ₂
Growth Time	Phase 1: 48 hours Phase 2: 4 to 6 days

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Method 59

Media	See note below for important directions. Phase 1: SP4 Glucose Broth Phase 2: SP4 Glucose Agar
Temperature	Phase 1: 35°C Phase 2: 35°C
Atmosphere	Phase 1: Aerobic Phase 2: CO ₂ (Candle Jar)
Growth Time	Phase 1: 7 to 28 days Phase 2: 5 to 15 days

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Method 60

Media	See note below for important directions. Phase 1: Mycoplasma Broth Phase 2: Mycoplasma Agar
Temperature	Phase 1: 35°C Phase 2: 35°C
Atmosphere	Phase 1: Aerobic Phase 2: 5 to 7% CO ₂
Growth Time	Phase 1: 48 hours Phase 2: 3 to 7 days

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Method 61

Media	See note below for important directions. Phase 1: SP4 Urea Broth Phase 2: A8 Agar
Temperature	Phase 1: 35°C Phase 2: 35°C
Atmosphere	Phase 1: Aerobic Phase 2: Anaerobic
Growth Time	Phase 1: 48 hours Phase 2: 4 to 6 days

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Method 65

Media	Nonselective Sheep Blood Agar or Columbia CNA with 5% Sheep Blood
Temperature	35°C
Atmosphere	5 to 7% Carbon Dioxide
Growth Time	48 hours

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Method 66

Media	Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	5% Carbon Dioxide
Growth Time	24 to 72 hours

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Method 67

Media	Tryptic Soy Agar (Soybean Casein Agar) or Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

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Method 68

Media	Chocolate Agar
Temperature	35°C
Atmosphere	Microaerophilic
Growth Time	48 hours

Method 69

Media	Malt Extract Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	5 to 7 days

Method 70

Media	Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

Method 71

Media	Potato Dextrose Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	7-14 days

Method 72

Media	Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 hours

Method 73

Media	Sabouraud Dextrose Emmons Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	2 to 4 weeks

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Method 74

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), Nonselective Sheep Blood Agar or Standard Methods Agar (Plate Count Agar)
Temperature	30°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

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Method 75

Media	Tryptic Soy Agar (Soybean Casein Agar) or Nonselective Sheep Blood Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

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Method 76

Media	Potato Dextrose Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	One week

Method 77

Media	Potato Dextrose Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	2 to 4 days

Method 78

Media	Sabouraud Dextrose Emmons Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	48 to 72 hours

Method 79

Media	Potato Dextrose Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	3 to 7 days

Method 80

Media	Nutrient
Temperature	25°C
Atmosphere	Aerobic
Growth Time	3 to 5 days

Method 81

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), Nonselective Sheep Blood Agar, Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	4 to 5 days

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